ENA CEM

(GB) ENGLISH Ena Cem instructions	pag. 2
(I) ITALIANO Ena Cem istruzioni	pag. 5
(D) DEUTSCH Ena Cem Gebrauchsanweisung	pag. 8
(F) FRANCAIS Ena Cem instructions	pag.13
(E) ESPAÑOL Ena Cem instrucciones	pag.17
(P) PORTUGUÊS Ena Cem instruções	pag.21
(CZ) ČESKY Ena Cem instrukce	pag.25
(RU) РУССКИЙ Ena Cem инструкции	pag.29

ENA
CEM

(GB) ENA CEM

Ena Cem is a dual curing radiopaque fluorescent luting composite, in dentine colour (UD2), for cementation of posts, ceramic and composite inlay, onlay veneers and crowns. It follows the standard ISO 4049:2000.

Composition

- MONOMER MATRIX: Diurethandimethacrylate, Iso-propyliden-bis (2(3)-hydroxy-3(2)-4(phenoxy)propyl)-bis(methacrylate)(Bis-GMA); 1,4 - Butandioldimethacrylate.
- TOTAL CONTENT OF FILLERS: 63% weight; glass filler: mean particle size 4,3 µm; highly dispersed silicone dioxide: mean particle size 0,012 µm

Indications: Ena Cem is used for cementation of:

- Post
- Ceramic laminated veneers, inlays, onlays, jacket crowns
- Composite inlays, onlays, crowns with thickness over 2 mm. (below 2 mm use Enamel plus HFO Hybrid light curing composite).

Contra-indications

Uncured resin could cause skin allergy: User should use gloves. In case of known allergy to some of the components do not use it.

Side effects

In deep cavities we suggest to use a liner in order to avoid pulpal reactions.

Materials to be avoided

Materials containing phenolics (like eugenol) and self etching primers could inhibit composite curing. Avoid the use of these materials as liners. This also applies to self-etching primers, which could prevent a proper curing of the cement. Therefore, the use of such primers should be foregone. We recommend the use of e. g. Ena Bond.

DIRECTION FOR USE

- | | |
|------------------------|----------|
| Working time: | 4 min. |
| Complete setting time: | 7-9 min. |

POST CEMENTATION

A. Complete the endodontic therapy, and prepare the tooth for the restoration by evaluating if a post is necessary, based on residual quantity of intact tooth. The loss of one or both interproximal walls in anteriors, of one or both wall mesial or distal in posteriors, suggest you to use a post, in order to better stabilize the restoration.

B. Make a radiograph to determine the appropriate diameter and depth of post space preparation. The post diameter should be slightly higher than the endodontic preparation. The post should be placed at same radicular depth as the height of dentinal body (for composite aesthetic restauration) or of post and core recommend the use of rubber dam.

C. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended.

D. Select the EnaPost Drill diameter which corresponds to the last Gates-Glidden drill to begin preparing the post space.

Use Drills with a slow-speed contra-angle with water spray in order to avoid tissues overheating. The drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove any tooth debris from canal by irrigating post space with water spray. Sequentially step up to the next larger EnaPost Drill until the preplanned diameter and depth are achieved.

E. Select the EnaPost that corresponds to the last drill used to prepare the post space and insert post into post space. For a better fitting shorten the post as necessary from either the apical or occlusal end (as clinical judgement dictates) with a Diamond separating disc, wetting the post with water to avoid fibers overheating.

F. Remove the post and gently clean it with an alcohol wipe and apply on the surface a mixture of Ena Bond with Ena Bond Catalyst (one drop+one drop).

G. Etch the cavity (coronal and radicular preparation) with Ena Etch 37% phosphoric acid for 2 minutes. For ideal adesion, before etching, a microblasting on cavity surfaces is recommended in order to clean and eliminate endodontic material debris.

H. Wash accurately the canal with a syringe to remove completely the acid. Suck water and dry the canal with paper point; to not dry with air in order to maintain wet the dentine.

I. Apply in the cavity and in the canal the mixture of Ena Bond and Ena Bond Catalyst. The adhesive should be wiped on the surface with a disposable microbrush or with a paper point. **Attention:** Ensure that microbrush reaches into the depths of the canal and that the Bonder is evenly rubbed in everywhere. Microbrush should not touch the surface or possibly get jammed. Dry with air to eliminate water and solvent residuals. Insert the post to check the canal and better push the adhesive in dentine tubules.

J. Apply a disposable automixing tip and extract Ena Cem dual composite cement from the syringe: the two components will be automatically mixed. Always discard from the syringe the first ca. 0.5 g of material. Place the cement directly into the canal filling it from the bottom of the cavity to the surface using a disposable endodontic applicator (Hawe-Neos); do not use any lentulo or rotating instrument. Apply some cement on the post surface and insert slowly the post to full depth, allowing excess cement to vent. Find the right position and, maintaining the post in position, wipe away any excess cement. **Note:** after each use remove the mixing tip and apply again the cap to the syringe.

K. Light cure for 60 seconds and proceed to the restoration using Enamel Plus HFO Mycro-Hybrid composite, both for aesthetic direct restoration and post and core, that could be immediately prepared once the restorative procedure has been completed.

Note: For further technical instructions please consult "EnaPost" manual.

LUTING OF PROSTHETIC APPLIANCE

Remove the temporary inlay and clean the cavity. Apply the rubber dam. Clean and dry the surface of the preparation. Try-in the appliance carefully. Etch tooth and apply the bonding. For ideal adhesion, before etching, a micro-blasting on cavity surfaces is recommended in order to clean and eliminate endodontic material debris.

Luting of ceramic appliances: For ceramic appliances, etch internal part of build-up by using Enamel plus ETCH Hydrofluoric Acid 9,6% for 60 seconds, wash with water, dry and apply Enamel plus ETCH Seal. Apply Ena Cem on the build-up to be luted, position it on the tooth. Remove composite excess and cure for at least 60 seconds from each side of the tooth. Check the occlusion, finish and polish with Enamel plus Shiny system, using burs, strips and diamond pastes.

Luting of composite appliances. For cementation of composite appliances sandblast the internal part, clean with alcohol and apply Ena Bond without curing it. For inlay, onlays and Veneers use Enamel plus HFO Hybrid light curing composite (we suggest OW or UD2 colour, follow Enamel plus HFO instructions). In case of inlay thickness over 2 mm use Ena Cem dual luting composite (see paragraph J and K).

Curing information

It's necessary to use a light-curing unit with a spectrum of 350 - 500 nm. The required physical results can be reached only if using multi walls reflecting unit. For this reason we suggest a periodical check of the light intensity following the manufacturer's instructions.

USE AND STORAGE

Store the product between 2°C and 10°C in a refrigerator. Do not use the product after the expiration date (see label on syringe). Due to hygienic reasons application tips should be used only once. The luting composite shows no chemical affinity to metals. Use the material at room temperature. Medical device, for dental use only: keep away from children. After each use remove the mixing tip and apply again the cap to the syringe. Avoid direct sunlight exposure.

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