# ENA Soft FLOW

## (GB) ENA Soft FLOW

ENA Soft FLOW is a light-curing temporary filling material. ENA Soft FLOW is an elastic material for temporary restorations. It has advantages of easy removal in one piece without damage to preparation. ENA Soft FLOW is suitable also for long term temporaries because it perfectly maintains the shape of tooth.

### ADVANTAGES

- simple to apply; not sticky;
- easy to remove in one piece;
- suitable for inlay & onlay;
- strong and long lasting fillings.

### Composition

- Diurethane dimethacrylate, mixture of isomers; 2-Hydroxyethyl methacrylate. Camphorquinone

### Indications

ENA Soft FLOW is used for:

- temporary fillings, especially for inlay & onlay

### Contra-indications

Uncured resin could cause skin allergy: user should use gloves. In case of known allergy to some of the components do not use it.

### Side effects

Inhalation: In case of inhalation, immediately consult with a medical specialist.

Skin contact: Wash off immediately with soap and enough water, and then consult with a dermatologist.

Eye contact: Rinse the eyes immediately with enough water, and then pour clean water into the eyelids more than 15 minutes, and then consult with an ophthalmologist.

### Materials to be avoided

Materials containing phenolics (like eugenol) could inhibit composite curing.

Avoid the use of these materials as liners.

### Instructions for use

- 1) Unscrew the cap and place a disposal application needle (minimum 1,2mm diam.);
- 2) Prepare the cavity.
- 3) Clean the cavity and dry it.
- 4) Apply sufficient amount of ENA Soft FLOW, depending on size of restoration.
- 5) Check occlusion and then remove excess material.
- 6) Light-cure for 10 seconds, for a 4 mm thickness, and for 20 seconds, for a 8 mm thickness, using a led (like CLEDPLUS – Micerium) or halogen light-curing unit with more than 1000 mW/cm<sup>2</sup> power (if between 550 mW/cm<sup>2</sup> and 1000 mW/cm<sup>2</sup> double the time).
- 7) The restoration can be removed with an instrument: it comes out in one piece.

### Use and storage

Do not store below 2°C and above 27°C. Do not use the product after the expiration date (see label on syringe).

Avoid package contamination. Use the material at room temperature.

Medical device, for dental use only: keep away from children. Turn back the spindle after taking out the material, to avoid sticking of the material, and close the container. After use, close container with cap and keep it closed.

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