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Additional Instructions

Ena Post is a POST made of glass fibers and resin used to reinforce restorations or core build-up in devitalized teeth.

It is available in 2 conicities (2% and 10%) and 4 sizes each.

They are used after preparing the root canal with a colour-coded drill calibrated with the post.

Indications

Ena Post is used for:

- Reinforcement of aesthetic restorations in devitalized teeth
- Reinforcement of core build-up in devitalized teeth

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REPROCESSING INSTRUCTIONS FOR REUSABLE MEDICAL DEVICES (ISO EN 17664)

Device(s): calibrated drills for Ena Post.

WARNINGS

Instruments are supplied not sterile and should be sterilized before the first use and after each use. Appropriate handling precautions should be taken during use and reprocessing. Gloves, eye protection and a mask should be worn.

Limitations on reprocessing

Repeated processing has minimal effect on these instruments; the working part is deteriorating mainly with use.

INSTRUCTIONS

Point of use: Remove excess soil under running water.

Container: No particular requirements. They should be protected from damage to the working part: burs holder is recommended

Transportation: It is recommended that instruments are reprocessed as soon as possible after using them. If transported wet there is an increased chance of staining or corrosion. Prolonged storage in disinfectant solutions may result in corrosion and should be avoided, however, these solutions must be validated for use with rotating instruments.

Cleaning:

- **automated cleaning:** Ultrasonic Cleaner (e.g. Biosonic) with general purpose liquid type UC30 or other liquids validated for cleaning of dental rotary instruments (manufacturers' instructions must be followed). Specific disinfectant solutions for ultrasonic can also be used only if validated for use with rotating instruments, specifically, for brushes and rubber points.

1. Load burs directly in a beaker or in a burs holder.

2. Run cycle, minimum 15 minutes.

3. When unloading check instruments and tips for complete removal of visible soil. If necessary repeat cycle or use manual cleaning.

- **manual cleaning:** 1. Rinse soil in excess from instrument under running cold water. 2. keep the instruments immersed in a sink reserved for the purpose in which case a disinfectant agent validated for processing of dental rotary instruments must be used and the disinfectant manufacturers' instructions must be followed and use wire brushes to remove any contamination. Care should be taken to avoid spreading contaminants by spraying or splashing during the brushing process 3. Rinse under clean running water or, even better, with ionized water.

Disinfection: Disinfectant solution may be used in accordance with label instructions. Only disinfectant agent validated for processing of dental rotary instruments must be used.

Drying: When drying is achieved as part of a washer disinfectant cycle do not exceed 140°C.

Maintenance: No particular requirements

Inspection and Function - Testing: After cleaning, inspect all instruments, with the aid of magnification if necessary, to ensure that all contamination has been removed. Repeat the cleaning process if necessary. Visually inspect for damage and eliminate broken or deteriorated instruments.

Packaging: single: a standard packaging material may be used. Instruments may be loaded into dedicated bur holder, wrap bur holder using appropriate method.

Sterilization: Vacuum autoclave (ISO17665). Cycle: 3 minutes minimum (18 minutes for diamond and carbide burs and Ena Post drills), 134°C maximum, 2,1 bar (follow autoclave manufacturers instructions).

Storage: No specific limitations related to storing duration. Storage should be in dry, clean conditions and at ambient temperature containers or pouches must be dry before opening to avoid recontamination of the contents with water.

Additional Information: When sterilizing multiple instruments with one autoclave cycle ensure that the sterilizer's maximum load, suggested by the manufacturer, is not exceeded.

Manufacturer contact: See brochure for telephone and address of local representative or telephone (+39) 0185 7887880

The instructions provided above have been validated by the medical device manufacturer as being CAPABLE of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the processing as actually performed using equipment, materials and personnel in the processing facility achieve the desired result. This requires validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

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